

USER INSTRUCTIONS

CAUTION

**SUNLAMP--DANGER--ULTRAVIOLET RADIATION
USE ONLY IN FIXTURE EQUIPPED WITH A TIMER**

**DANGER--ULTRAVIOLET RADIATION
FOLLOW INSTRUCTIONS**

Avoid overexposure. As with natural sunlight, overexposure can cause eye and skin injury and allergic reactions. Repeated exposure may cause premature aging of the skin and skin cancer. **WEAR PROTECTIVE EYEWEAR; FAILURE TO MAY RESULT IN SEVERE BURNS OR LONG-TERM INJURY TO THE EYES.**

Medications or cosmetics may increase your sensitivity to the ultraviolet radiation. Consult physician before using sunlamp if you are using medications or have a history of skin problems or believe yourself especially sensitive to sunlight. If you do not tan in the sun, you are unlikely to tan from the use of this product.

The instructions accompanying the sunlamp product should ALWAYS be followed to avoid potential injury.

ALWAYS RETAIN THIS FORM FOR FUTURE REFERENCE AND LAMP COMPATIBILITY INFORMATION.

<u>MODEL</u>	<u>COMPATIBLE WITH</u>	<u>LENGTH</u>	<u>WATTS</u>	<u>BASE</u>
SS71/BP Super Sun® II	-- Wolff Bellarium S SA1-71-100W Wolff System Diamond Sun S DSS71-T12-100W Wolff System Speed 175 SPD 75-T12 Wolff System Speed 205 SPD 71-T12-205 Wolff System Accelerator ACC71-T12-100W Preheat Bipin Wolff System EverGlo 71"/100W Wolff Crystal Sun S DA1-763-100W BiPin Wolff Velocity VEL 71-T12-100W Wolff System Body Scan X Plus; F71/100W Bipin; Preheat ERS Sol-Soft AH3-24-100W Tan America VIP Gold F71/2.6/HO/BP, Part 46132	71"	100	Bi-Pin
SS71/BP 2.6 Super Sun®	-- Voltarc F71T12/BL/BP 2.6 Wolff Crystal Sun S DA1-763-100W Bi-Pin Sylvania Cosmolux S CA1-12-100W Cosmedico Cosmostar Longlife 12503/100W Preheat-BiPin	71"	100	Bi-Pin
SS71/VHO 2.6 Super Sun®	-- Sunfit Professional VX 160W	71"	160	Bi-Pin
SS71/BP 3.5 Super Sun®	-- Puretán F71T12-2.6-BP Wolff Crystal Sun S DA1-763-100W BiPin Cosmedico Goldarium S SA1-12-100W ERS Sol-Soft AH3-24-100W Philips Cleo Swift F71T12 TL 100W Sunfit Professional X 100W Tan America VIP Gold F71/2.6/HO/BP, Part 46132	71"	100	Bi-Pin

INTERLECTRIC CORPORATION
1401 LEXINGTON AVENUE, WARREN, PA 16365



9/16/10

Super Sun® Lamps
1-800-722-2184 (PA: 814-723-6061)

CONTINUED ON OTHER SIDE

USER INSTRUCTIONS

CONTINUED FROM OTHER SIDE

<u>MODEL</u>	<u>COMPATIBLE WITH</u>	<u>LENGTH</u>	<u>WATTS</u>	<u>BASE</u>
SS71/BP 5.5 Super Sun®	-- Light Sources F71T12-BL-BP/42 Puretán F71T12-4.2 BPS+ Voltarc F71T12-BL 42 PINK-BP Voltarc Euro Bronze Magnifier 70+; F71T12-BiPin+ Voltarc Genesis 105W Plus F71T12-BP-105W Sunbronze F71T12-BL-BP/44 Wolff System EverGlo 71" 100W Wolff System Speed 215 SPD71-T12-215 Wolff System Speed 205; SPD71-T12-205; BiPin Wolff System Diamond Sun S DSS71-T12-100W Preheat Bipin Wolff System, Velocity Extreme; VE71-T12-100W BiPin Wolff System Velocity Plus VLP71-T12-100/120W Bi-Pin Wolff System Velocity Select VS71-T12-100W Bi-Pin Intensity IT 71-T12-100W Bi-Pin Interlectric Super Sun SS71/BP 4.4	71"	100	Bi-Pin
SR71/BP Plus Super Sun® Reflector	-- Philips R-UVA 100W Philips TL 100W/09R R-UVA Sunfit RS RUVA-S-Highpower 100W	71"	100	Bi-Pin
SR71/BP/GEM-R Super Sun®	-- Sunfit RX 100W Professional Philips TL 100W Professional Cosmedico, Goldarium SR SA1-15-100W	71"	100	Bi-Pin
SS71/BP/A-RAY Super Sun®	-- Philips TL 100W/09N UV-A Sunfit S Super-12 Highpower 100W	71"	100	Bi-Pin
SS71/BP/3.2 Super Sun®	-- JK Soltron S 100W Cosmolux AS A1-12-100W Preheat Bi-Pin Wolff System Diamond Sun DS71-T12-100W Preheat Bi-Pin	71"	100	Bi-Pin
SS71/BP/PANTHER® Super Sun®	-- Cosmedico Hi-Tan Preheat Bi-Pin 100W Radiance 7000 F71	71"	100	Bi-Pin
SS71/BP/PHOENIX® Super Sun®	-- Sunbronze F71T12-BL-BP/85 EuroSun S3 ES3-71-100W Cosmedico Hi-Tan Preheat Bi-Pin 100W Radiance 7000 F71	71"	100	Bi-Pin
SS71/BP/PHOENIX® Extreme Super Sun®	-- Cosmedico Hi-Tan Preheat Bi-Pin 100W Radiance 7000 F71	71"	100	Bi-Pin

INTERLECTRIC CORPORATION
1401 LEXINGTON AVENUE, WARREN, PA 16365



9/16/10

Super Sun® Lamps
1-800-722-2184 (PA: 814-723-6061)